



***Scintillation Materials, Detectors and Electronics***

## **OPERATING MANUAL**

### ***DETECTORS FOR MEDICAL IMAGING SYSTEMS***

#### **I. DESCRIPTION**

Detector for medical imaging systems based on sodium iodide (NaI) polycrystal (hereinafter referred to as "detector") is used for registration and gamma radiation spectrometry in the photons energy range from 60 keV 360 keV for medical imaging systems.

Detector consists of scintillation crystal and protective glass encapsulated in aluminum housing. Scintillation crystal is wrapped in reflector material for maximum light detection.

#### **II. UNPACKING AND PRODUCT INSPECTION**

If on receipt of the products external damages of the packing are discovered, one should confirm the damages with the assistance of the shipping agent. The manufacturer can not replace the damaged products without documentarily confirmed damage of the packing.

It is allowed to unseal the packing only after isolation under normal climatic conditions no less than 24 hours.

Take out the certificate (certificate of technical characteristics inspection of the detector). Take out the detector and additional parts (if any), and set it in a safe place.

Make sure that all elements specified in the packing list are included into the packing. If there are several detectors in the packing, it is necessary to define the correspondence of certificates by serial numbers to each of the detectors.

The buyer shall carry out primary visual control of products in accordance with the specification within 30 (thirty) days from the date of products receipt as follows:

1. Check of appearance (any mechanical damages on the output window and housing of the container, indentations and corrosion on the metallic surface of the container; facets and foliations on the glue line, hydrolytic film sections of white, yellow-green color that is considered to be scintillator destruction resulting its chemical interaction with moisture; ungluing of output window from scintillator for more than 4mm<sup>2</sup>; chips and cracks in scintillator with dimensions more than 3 mm are not acceptable).

Appearance of detectors is controlled by visual inspection without using of magnifying facilities at the light from 10 lux to 50 lux.

Chips and cracks of smaller dimensions are not determined. Cracks and foliations of reflective envelope are considered to be construction integrity damaging.

2. Verification of the detectors marking (information about the type of the detectors, serial number of the detectors, date of production).



If at the primary visual control any of the above mentioned visual defects is not found in products, it is considered as having passed visual control and is suitable for testing.

In case of discovery of visual defects or other disparities the customer should turn in written form for explanations to the QC department of the manufacturer (QC@isc.kharkov.com) within 30 (thirty) days from the date of products receipt.

If after expiration of 30 (thirty) days from the date of products receipt the customer does not provide any information about the discovered visual defects during the primary visual control the manufacturer should not accept claims on disparity of the appearance of the products.

In case irreparable defects (mechanical damages on the output window and housing of the container, indentations and corrosion on the metallic surface of the container; facets and foliations on the glue line, hydrolytic film sections of white, yellow-green color that is considered to be scintillator destruction resulting from its chemical interaction with moisture; ungluing of output window from scintillator for more than 30 %; chips and cracks in scintillator with dimensions more than 3 mm.) are found in the products, they should be returned to the manufacturer for warranty repair or replacement. It is necessary to provide packing with the use of original packing materials by which it is possible to make certain that no damage occurred during transportation and storage. It is also necessary to attach the proper warning labels on the packing to provide careful handling of the products.

The products should be returned with the following accompanying documents:

- brief description of reasons for return;
- name of the buyer and contact telephone number.

### III. TECHNICAL SPECIFICATION

Detector for medical imaging systems is produced according to the design documentation.

Crystal material .....	NaI (TI)
Density, g/cm <sup>3</sup> .....	3.67
Hygroscopic .....	high
Energy range .....	60 – 360 keV
Energy resolution for 122 keV ( <sup>57</sup> Co), % with 48 mm PMT photocathode diameter	
with the thickness of crystal 6.5 mm .....	≤ 14
with the thickness of crystal 9.5 mm .....	≤ 15
Energy resolution for 59.5 keV ( <sup>241</sup> Am), % with 48 mm PMT photocathode diameter	
with the thickness of crystal 6.5 mm .....	≤ 17
with the thickness of crystal 9.5 mm .....	≤ 19
Operating temperature range, °C .....	+15 ...+30
Physical load .....	no more than 3 Hz with acceleration 3g



#### **IV. TESTING PROCEDURE**

1. The measurements should be carried out in enclosed camera under normal climatic conditions:
  - air temperature from 15 to 35°C;
  - relative air humidity from 45 to 80%;
  - atmosphere pressure from 84 to 106 kPa (from 630 to 800 mm Ht).
2. For registration of pulse height spectrum pulse height analyzers with the number of channels 1024 or 4096 with output recording device of any type which allows representing pulse height spectrum in the form convenient for processing should be used.

All equipment should be grounded.

The measurements should be carried out not earlier than for 30 minutes after the last device is enabled.

3. As gamma-radiation source sealed radioisotopic source of ionizing radiation with certain energy is used (basically <sup>57</sup>Co or <sup>241</sup>Am).

For high-energy radiation sources intrinsic width of ionizing radiation source line should correspond to indicated in the specification and should be smaller than pulse height distribution curve for one-two degrees.

For measurements of detectors for gamma-cameras radiation source with one-hole collimator is used.

4. Before measurements are made detector is stored in darkness within 12 hours.
5. Before measurements are made PMT is stored in darkness under high voltage within the time PMT is stabilized under operating conditions.

All characteristics are measured in a completely dark state of the scintillator and PMT.

6. Detector is installed with output window upwards. Then PMT is installed according to the scheme specified in quality document (the number of points and their distribution can be different as agreed by the customer).
7. Opposite the output window of detector collimated gamma-radiation source is installed in such way that PMT and source axes are overlapped. Optical contact between detector and PMT is provided by paraffinic oil (or another immersion oil) which is cleaned with ethyl alcohol from the output window of the detector and PMT after measurements are made.

At spectrum collection integral pulse count rate shouldn't exceed 3,0 thousand counts per second reached by fitting of the distance between the source of ionizing radiation and detector or fitting of source activity.

8. After high voltage supply and establishment of operating conditions for PMT pulse height spectrum is taken. Pulse height spectrum is recorded within the time not more than 3,0 thousand pulses are collected in the channel corresponding to full absorption peak maximum.
9. Processing of the measurements results. The amplitude resolution is evaluated according to the following equation:

$$R = \frac{\Delta V}{V} * 100,$$

where:  $\Delta V$  – the width of full absorption peak at FWHM;

$V$  – amplitude of the signal.

The detector passed the testing if measurement data are in the range of  $\pm 5\%$  of the results for the last test at the manufacturer's company indicated in the certificate.



## **V. SAFETY CONDITIONS AND OPERATING RULES**

The product is not toxic, explosive or radioactive.

Conditions for normal operation:

1. Do not drop the product.
2. For indoor use or out door use (in dry medium).
3. Relative air humidity should not exceed 95%.
4. Physical load should not exceed 50 Hz.
5. Temperature range is from -20°C to 40°C.

In case of damage or failure the in-house repair of the product is not allowed. Damage product should be hermetically packed in polyethylene and sent to the manufacturer.

The product should be cleaned outside with a damp cloth and ethyl alcohol. The product should be cleaned only when it is out of running. Do not dip the product into the water.

## **VI. TRANSPORTATION AND STORAGE**

The transportation of the products in the manufacturer's packing can be carried out by all means of transport to any distance with the observance of the following conditions:

1. The transportation should be carried out in clean covered places.
2. The installation and the fastening of the boxes in the transport should ensure their stable position eliminating the possibility of displacement and collision between themselves and the sides of the transport.
3. Warning marking should be observed on every stage of the transportation.
4. Avoid the fall of the boxes with the products.

The products in the manufacturer's packing should be stored indoors with good ventilation at the temperature from 15°C to 35°C (from 60°F to 95°F) and relative air humidity not exceeding 80% at the absence of corrosive impurities in the environment.

The products should be protected from huge temperature swings.

## **VIII. WARRANTY AND MAINTENANCE SERVICE**

The manufacturer warrants that the products should be free of defects in materials used for manufacturing or production process.

The manufacturer warrants the conformity to the products quality at the observance by the customer of transportation and storage conditions, and the operation rules established by this operation manual:

- 1) Storage warranty period is 30 (thirty) months from the date of the final quality control of the products at the manufacturer's enterprise.
- 2) Operation warranty period is 10 000 (ten thousand) hours from the date of the products putting into operation within the storage warranty period.

The warranty applies only to the original customers of the product and only to the product with serial number still legible.



**Scintillation Materials, Detectors  
& Electronics**

The warranty does not apply to the product whose performance fails as a result of misuse, mishandling, abuse, accident, physical damage or improper installation, or submersion in water. The warranty is not valid if the product has been modified or altered without manufacturer's express approval.

The manufacturer's obligation shall be limited to the repair or replacement of the product, at manufacturer's option, which has been rejected as defective within the warranty period.

The manufacturer is not responsible for damages including incidental or consequential damages resulting from improper use or failure of products operation under this warranty.

Any claim on the warranty should be made by the original customer within the warranty period. The customer should apply to the Quality Control Department (QC@isc.kharkov.com) in order to receive additional information.

All products returning for warranty repair or replacement should be packed with the use of original packing materials by which it is possible to make certain that no damage occurred during transportation and storage.

All products returning for warranty repair or replacement should be accompanied by the documents confirming the purchase, as well as by rejection sheet, which confirms that use of the product is not possible through the manufacturer's fault.

The products replaced under warranty repair, are not subject to the payment

Returned products are found out to be in good operational state or damaged for causes which are not covered by this warranty will be subject to a service charge. After the repair is finished the products will be returned to the customer freight prepaid.